

Anti-Smoothelin antibody [R4A] (STJ16100219)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Dog/Chicken/Feline/Monkey/Swine

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	6525
Gene Symbol	SMTN
UniProt ID	SMTN_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

