

Anti-NSP-A and NSP-B antibody [RNL-2] (STJ16100220)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ICC/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Monkey/Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:25-1:100 ICC IHC-F 1:50 IHC-p 1:50 WB 1:100-1:1000
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

