

Anti-NSP-A and NSP-B antibody [RNL-3] (STJ16100221)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ICC/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Monkey/Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:100-1:200 WB 1:100-1:500 ICC 1:100-1:200 IHC 1:100-1:200
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

