

## Anti-Basal cell Cytokeratin antibody [RCK103] (STJ16100229)

### GENERAL INFORMATION

|               |   |
|---------------|---|
| Product Type  | Primary antibodies  |
| Applications  | FC/ICC/IHC-F/WB   |
| Host / Source | Mouse   |
| Reactivity    | Human/Dog/Chicken/Guinea Pig/Hamster/Quail/Rabbit/Rat/Swine/Zebrafish |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1 mg/mL  |
| Conjugation         | Unconjugated   |
| Dilution Range      | FC 1:00-1:200<br>ICC<br>IHC-F 1:100-1:200<br>WB 1:100-1:1000 |
| Formulation         | PBS with 0.09% sodium azide                                  |
| Isotype             | IgG1   |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.               |

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

