

## Anti-Pan-Keratin antibody [80] (STJ16100232)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-F/WB/FC        |
| Host / Source | Mouse              |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |                                                |
|---------------------|------------------------------------------------|
| Clonality           | Monoclonal                                     |
| Concentration       | 15 Mu g/mL                                     |
| Conjugation         | Unconjugated                                   |
| Dilution Range      | IHC 1:50<br>IHC 1:10-1:50                      |
| Formulation         | Culture medium with 0.01M Sodium Azide         |
| Isotype             | IgG1                                           |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

