

Anti-Vimentin antibody [V9] (STJ16100234)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F/WB/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	25 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:300 IHC-P 1:50-1:100
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7431
Gene Symbol	VIM
UniProt ID	VIME_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

