

Anti-Glial Fibrillary Acidic Protein antibody [GF-01] (STJ16100238)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IP/WB/ICC
Host / Source	Mouse
Reactivity	Pig/Cat/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	2670
Gene Symbol	GFAP
UniProt ID	GFAP_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

