

Anti-Keratin 4 antibody [6B10] (STJ16100240)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	15 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:40 IHC-P 1:10
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	3851
Gene Symbol	KRT4
UniProt ID	K2C4_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

