

## Anti-Keratin 8 antibody [M20] (STJ16100241)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB/FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	15 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:20 WB FC
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3856</a>
Gene Symbol	<a href="#">KRT8</a>
UniProt ID	<a href="#">K2C8_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

