

Anti-Keratin 13 antibody [2D7] (STJ16100243)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	15 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-P 1:5-1:25
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	3860
Gene Symbol	KRT13
UniProt ID	K1C13_HUMAN

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

