

## Anti-Keratin 10 antibody [DE-K10] (STJ16100249)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	WB 1:100-1:250 IHC 1:5-1:20
Formulation	Contains 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3858</a>
Gene Symbol	<a href="#">KRT10</a>
UniProt ID	<a href="#">K1C10_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---