

Anti-Neurofilament antibody [403] (STJ16100255)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/IF/WB
Host / Source	Mouse
Reactivity	Bovine/Guinea Pig/Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	ICC 1:10 IF 1:10 IHC 1:10
Formulation	Contains 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------