

## Anti-Keratin 14 antibody [RCK107] (STJ16100264)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Human/Dog/Rat/Swine

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:100-1:200 ICC IHC-F 1:100-1:200 WB 1:100-1:1000
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3861</a>
Gene Symbol	<a href="#">KRT14</a>
UniProt ID	<a href="#">K1C14_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

