

Anti-TLR3 antibody [TLR3.7] (STJ16100285)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/IF/IP/IHC-P
Host / Source	Mouse
Reactivity	Human/Dog/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 µg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:50 IP 1:50 IHC 1:50 IP 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
UniProt ID	TLR3_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---