

Anti-S100A8/A9 antibody [2+70E+11] (STJ16100301)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/ELISA/IF/IP/WB
Host / Source	Mouse
Reactivity	Human/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	FC 1:50 WB 1:50 IHC 1:50 IP 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

