

Anti-Macrophages antibody [25F9] (STJ16100303)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	FC 1:10 WB 1:10 IHC 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---