

## Anti-IL2RA antibody [4C9] (STJ16100332)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Dilution Range      | IHC 1:10 – 1:50   |
| Formulation         | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Isotype             | IgG2b   |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.                    |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">3559</a>        |
| Gene Symbol | <a href="#">IL2RA</a>       |
| UniProt ID  | <a href="#">IL2RA_HUMAN</a> |

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

