

## Anti-EMA antibody [E29] (STJ16100362)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                   |
|---------------------|-------------------------------------------------------------------|
| Clonality           | Monoclonal                                                        |
| Conjugation         | Unconjugated                                                      |
| Dilution Range      | IHC 1:100 – 1:500                                                 |
| Formulation         | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Isotype             | IgG2ak/IgG1k                                                      |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.                    |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">4582</a>       |
| Gene Symbol | <a href="#">MUC1</a>       |
| UniProt ID  | <a href="#">MUC1_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

