

## Anti-MadCAM-1 antibody [314G8] (STJ16100465)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA/IHC-P/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Dilution Range	WB 1:50 FC 1:50 IHC 1:50 ELISA 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">8174</a>
Gene Symbol	<a href="#">MADCAM1</a>
UniProt ID	<a href="#">MADCA_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---