

## Anti-Myosin antibody [414] (STJ16100469)

### GENERAL INFORMATION

|               |                            |
|---------------|----------------------------|
| Product Type  | Primary antibodies         |
| Applications  | IHC-F/IHC-P/IF/WB          |
| Host / Source | Mouse                      |
| Reactivity    | Chicken/Human/Mouse/Rabbit |

### PRODUCT PROPERTIES

|                     |                                                |
|---------------------|------------------------------------------------|
| Clonality           | Monoclonal                                     |
| Conjugation         | Unconjugated                                   |
| Dilution Range      | IF 1:10<br>IHC 1:5-1:10                        |
| Formulation         | Contains 0.09% sodium azide                    |
| Isotype             | IgM                                            |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|