

## Anti-a-Tubulin antibody [TU-16] (STJ16100475)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP/WB/IHC-P/ICC/ELISA
Host / Source	Mouse
Reactivity	Mouse/Rat/Pig/Nicotiana/Human/Dog

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Tris buffered saline (TBS) solution with 15 mM sodium azide
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

