

Anti-Tubulin Alpha antibody [TU-02] (STJ16100478)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC
Host / Source	Mouse
Reactivity	Human/Mouse/Pig/PlantS./Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: not known; tracer: 0, 001-100 Mu g/ml for 30 min at RT). Immunohistology (formalin-fixed: 1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues is enhanced after boiling tissue sections in 10mM citrate buffer, pH 6.0, for
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgMk
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	7846
Gene Symbol	TUBA1A
UniProt ID	TBA1A_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

