

Anti-TLR4 antibody [HTA125] (STJ16100496)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ELISA/IF/IP/WB
Host / Source	Mouse
Reactivity	Human/Dog/Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7099
Gene Symbol	TLR4
UniProt ID	TLR4_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---