

## Anti-CD62L antibody [FMC46] {PE} (STJ16100507)

### GENERAL INFORMATION

|               |                                           |
|---------------|-------------------------------------------|
| Product Type  | Primary antibodies                        |
| Applications  | FC                                        |
| Host / Source | Mouse                                     |
| Reactivity    | Human/Rhesus/Bovine/Cynomolgus Monkey/Dog |

### PRODUCT PROPERTIES

|                     |                                                    |
|---------------------|----------------------------------------------------|
| Clonality           | Monoclonal                                         |
| Conjugation         | PE                                                 |
| Formulation         | PBS with 0.09% sodium azide and 1% BSA, 5% sucrose |
| Isotype             | IgG2b                                              |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.     |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

