

Anti-CD227 antibody [C595] (STJ16100530)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/FC/IHC-F/ELISA/IF
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.09% sodium azide
Isotype	IgG3
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

