

## Anti-CD226 antibody [DX11] {PE} (STJ16100571)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	PE
Formulation	PBS with 0.09% sodium azide and 1% BSA, 5% sucrose
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">10666</a>
Gene Symbol	<a href="#">CD226</a>
UniProt ID	<a href="#">CD226_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---