

Anti-Transferrin antibody [HTF-14] (STJ16100593)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ICC/ELISA
Host / Source	Mouse
Reactivity	Human/Pig/Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by protein-A affinity chromatography.
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7018
Gene Symbol	TF
UniProt ID	TRFE_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

