

Anti-Interleukin 6 antibody [8F4] (STJ16100598)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FUNC/ELISA/IP/WB
Host / Source	Mouse
Reactivity	Human/Cynomolgus Monkey/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 ug/mL
Conjugation	Unconjugated
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3569
Gene Symbol	IL6
UniProt ID	IL6_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---