

## Anti-Alpha-Fetoprotein antibody [AFP-11] (STJ16100609)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P/ELISA        |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 1 mg/mL  |
| Conjugation         | Unconjugated   |
| Purification        | Purified by protein-A affinity chromatography.                   |
| Formulation         | Phosphate buffered saline (PBS) solution with 15 mM sodium azide |
| Isotype             | IgG1   |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.                   |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">174</a>        |
| Gene Symbol | <a href="#">AFP</a>        |
| UniProt ID  | <a href="#">FETA_HUMAN</a> |

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

