

Mouse Anti-Horseradish/Tag HRP antibody {HRP} (STJ16100611)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	WB/ICC/ELISA
Host / Source	Mouse
Reactivity	Horseradish/Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	HRP
Dilution Range	WB ICC ELISA
Formulation	Phosphate buffered saline (PBS) solution with 15 mM Sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

