

Anti-CD90 antibody [F15-42-1] {FITC} (STJ16100674)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Cynomolgus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.1 mg/mL
Conjugation	FITC
Formulation	PBS with 0.09% sodium azide and 1% BSA
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---