

Anti-Liver Fatty Acid Binding Protein antibody [L2B10] (STJ16100686)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Baboon/Dog/Rat/Swine

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---