

Anti-Fyn antibody [FYN-1] (STJ16100691)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP/WB/IHC-P/ICC
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	2534
Gene Symbol	FYN
UniProt ID	FYN_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

