

## Anti-nuclear antigen 1 of Epstein-Barr virus antibody [E1-2.5] (STJ16100693)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/ELISA
Host / Source	Mouse
Reactivity	Epstein-Barr Virus

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---