

## Anti-MUC1/Polymorphic epithelial mucin antibody [2G7] (STJ16100705)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-F/IHC-P/IF/FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">4582</a>
Gene Symbol	<a href="#">MUC1</a>
UniProt ID	<a href="#">MUC1_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

