

Anti-MUC5AC/M1 antibody [9-13M1] (STJ16100711)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Mouse
Reactivity	Cat/Cow/Human/Monkey/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Immunoblotting (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues requires 0, 1% pronase digestion (0, 1% CaCl ₂ , 0, 1N NaOH) for
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4586
Gene Symbol	MUC5AC
UniProt ID	MUC5A_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

