

## Anti-MMP-9 antibody [2C3] (STJ16100713)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">4318</a>
Gene Symbol	<a href="#">MMP9</a>
UniProt ID	<a href="#">MMP9_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---