

## Anti-T-cell Receptor & Epithelial Tumours antibody [R73xCC52] (STJ16100721)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FUNC
Host / Source	Mouse
Reactivity	Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---