

Anti-Tumor Antigens of Epithelial Origin antibody [CC52] (STJ16100726)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F/FC/WB/IP/ELISA
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Culture supernatant with 0.05% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---