

Anti-CD21 antibody [4C5] (STJ16100747)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IF/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	60 µg/mL
Conjugation	Unconjugated
Dilution Range	IHC-F 1:20-1:50 IF FC 1:50
Formulation	serum with 0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---