

Anti-Granzyme B antibody [GrB-7] (STJ16100752)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/ELISA/IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Purification	Protein G purified
Dilution Range	ICC 1:20
Formulation	Serum free tissue culture supernatant with 0.1% BSA and 0.09% sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	3002
Gene Symbol	GZMB
UniProt ID	GRAB_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---