

## Anti-p21/WAF1 antibody [WA-1] (STJ16100762)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF/IHC
Host / Source	Mouse
Reactivity	Human/Chimpanzee/Monkey/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Flow Cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (1-2 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; requires boiling tissue sections in 10
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">1026</a>
Gene Symbol	<a href="#">CDKN1A</a>
UniProt ID	<a href="#">CDN1A_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

