

## Anti-Rabies Virus antibody [RAB-50] (STJ16100770)

### GENERAL INFORMATION

|               |                      |
|---------------|----------------------|
| Product Type  | Primary antibodies   |
| Applications  | ELISA/WB/IF/IHC/NEUT |
| Host / Source | Mouse                |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                                                                  |
| Concentration       | Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.                                                                                                                                                                                                          |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                                |
| Purification        | Affinity purified from tissue culture.                                                                                                                                                                                                                      |
| Dilution Range      | ELISA (solid phase: not known; tracer: 0.001-100 µg/ml for 30 min at RT). Immunoblotting (1.0-2.0 µg/ml). Immunofluorescence (1-2 µg/ml). Immunohistology (1-2 µg/ml for 330 min at RT; an appropriate antigen retrieval method for staining of formalin-fi |
| Formulation         | PBS with 0.02% Sodium Azide.                                                                                                                                                                                                                                |
| Isotype             | IgG2b Kappa                                                                                                                                                                                                                                                 |
| Storage Instruction | Store for up to 1 year at 2-8°C upon receipt.                                                                                                                                                                                                               |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

