

## Anti-LTA antibody [55] (STJ16100772)

### GENERAL INFORMATION

|               |                      |
|---------------|----------------------|
| Product Type  | Primary antibodies   |
| Applications  | IHC-F/FC/ELISA/IF/WB |
| Host / Source | Mouse                |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal                                     |
| Concentration       | > 200 Mu g/mL                                  |
| Conjugation         | Unconjugated                                   |
| Formulation         | PBS with 0.1% BSA and 0.02% sodium azide       |
| Isotype             | IgG3   |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">4049</a>       |
| Gene Symbol | <a href="#">LTA</a>        |
| UniProt ID  | <a href="#">TNFB_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|