

## Anti-C9 antibody (C-Term) [X197] (STJ16100776)

### GENERAL INFORMATION

|               |                         |
|---------------|-------------------------|
| Product Type  | Primary antibodies      |
| Applications  | IHC-F/FC/ELISA/IHC-P/WB |
| Host / Source | Mouse                   |
| Reactivity    | Human                   |

### PRODUCT PROPERTIES

|                     |                                                |
|---------------------|------------------------------------------------|
| Clonality           | Monoclonal                                     |
| Concentration       | 100 Mu g/mL                                    |
| Conjugation         | Unconjugated                                   |
| Formulation         | PBS with 0.1% BSA and 0.02% sodium azide       |
| Isotype             | IgG1                                           |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### TARGET INFORMATION

|                  |                           |
|------------------|---------------------------|
| Gene ID          | <a href="#">735</a>       |
| Gene Symbol      | <a href="#">C9</a>        |
| UniProt ID       | <a href="#">CO9_HUMAN</a> |
| Immunogen Region | C-Term                    |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|