

Anti-Retinol Binding Protein antibody [G4E4] (STJ16100778)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IHC
Host / Source	Mouse
Reactivity	Chimpanzee/Goat/Human/Monkey/Mouse/Rabbit/Rat/Rhesus

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: not known; tracer: 0, 001-100 µg/ml for 30 min at RT). Immunoblotting (0, 5-1, 0 µg/ml). Immunohistology (formalin-fixed: 2-4 µg/ml for 30 min at RT; staining of formalin-fixed tissues requires boiling tissue sections in 10mM citr
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	5950
Gene Symbol	RBP4
UniProt ID	RET4_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

