

Anti-PAI-1 antibody [MA-33H1F7] (STJ16100781)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FUNC/ELISA/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rabbit/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	5054
Gene Symbol	SERPINE1
UniProt ID	PAI1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---