

Anti-Tenascin-C antibody [EB2] (STJ16100797)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB/ELISA
Host / Source	Mouse
Reactivity	Human/Rabbit/Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 1% BSA & 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	3371
Gene Symbol	TNC
UniProt ID	TENA_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---