

Anti-Epidermal Growth Factor antibody [E5] (STJ16100805)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F/WB/ELISA
Host / Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu g/mL
Conjugation	Unconjugated
Dilution Range	IHC-P 1:5 IHC-F 1:5 WB ELISA
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---