

Anti-Bromodeoxyuridine antibody [IIB5] (STJ16100807)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC/IHC-F/IHC-P/ICC
Host / Source	Mouse
Reactivity	Other

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:100-1:200 IHC IHC-F IHC-P
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

